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**D 4010**

**Type 4739**

1. Insert batteries as shown. Uses two 1.5 V type AA (LR 6) batteries. For best performance, we recommend Duracell® Ultra alkaline batteries. Batteries will last approx. 3 months when brushing twice a day for two minutes.
2. Remove protection cap, wet brushhead and apply toothpaste.
3. Place brushhead on teeth and switch on.
4. Brush for at least 2 minutes (recommended by dentists). Guide the brushhead from tooth to tooth, first the outsides (a), then the insides (b), finally the chewing surfaces (c). During the first days of use, gums may bleed slightly. Should bleeding persist after 2 weeks, contact your dentist.
5. Replace brushhead when indicated (a = new, b/c = replace).
6. Clean as shown under running water, then dry.
7. Dispose of used batteries according to local waste collection requirements.

Adult supervision recommended when used by children.

Limited Warranty: In case of defects in materials or workmanship, the appliance will be replaced within 2 years of date of purchase. If the appliance exhibits such a defect call 1-800-272-8611 (US residents) or 1-800-387-6657 (Canadian residents) for address to return appliance. We will send you a replacement appliance free of charge. Proof of date of original purchase is required. This warranty does not cover use and wear of batteries and brushhead, nor does it cover damage to finishes, normal wear, or damage resulting from accident, misuse, dirt, tampering, or unreasonable use.

ALL IMPLIED WARRANTIES, INCLUDING ANY IMPLIED WARRANTY OF FITNESS FOR ANY PARTICULAR PURPOSE, ARE EXPRESSLY DISCLAIMED. IN NO EVENT WILL MANUFACTURER BE RESPONSIBLE FOR CONSEQUENTIAL, INCIDENTAL OR SPECIAL DAMAGES RESULTING FROM THE USE OF THIS PRODUCT.

Some states/provinces do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation may not apply to you.

This warranty gives you specific legal rights, and you may have other rights which vary from state to state/province to province.

**BRAUN**

is a trademark of Braun GmbH, Kronberg/Ts., Germany

**Oral-B**

and INDICATOR are trademarks of Oral-B Laboratories

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USA  
Distr. by:  
Gillette Commercial Operations  
North America  
One Gillette Park,  
Boston, MA 02127-1096

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Offer valid for the following Braun Oral-B rechargeable power toothbrushes only: 3D Excel, 3D, Plak Control Ultra, Plak Control & Plaque Remover. Offer not valid on Braun Oral-B Battery Toothbrush and brushheads.

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**To the dealer:**  
Upon receipt of this coupon toward the purchase of the specified product, we will reimburse you the face value of the coupon plus regular handling. Application for redemption on any other basis may constitute fraud and will, at our option, void coupon presented. Applications for reimbursement accepted from principals only.  
Mail to: Braun Oral-B, Box 3000, Saint John, N.B. E2L 4L3  
A reduction in any applicable taxes payable is included in the coupon face value.

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